

Product Datasheet

B9D2 antibody (orb1151124)

Description	Rabbit polyclonal antibody to B9D2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human B9D2(NP_085055.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19kDa
Uniprot ID	Q9BPU9
Dilution Range	IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.