

Product Datasheet

DBI antibody (orb1151125)

Description	Rabbit polyclonal antibody to DBI
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human DBI(NP_001171513.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11kDa
Uniprot ID	P07108
Dilution Range	WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.