

Product Datasheet

FYB antibody (orb1151126)

Description	Rabbit polyclonal antibody to FYB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human FYB(NP_001456.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120KDa
Uniprot ID	O15117
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.