

## Product Datasheet

### CENPA (Phospho-S7) antibody (orb1151129)

|                            |                                                                               |
|----------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CENPA (Phospho-S7)                              |
| <b>Species/Host</b>        | Rabbit                                                                        |
| <b>Reactivity</b>          | Human                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                  |
| <b>Tested Applications</b> | WB                                                                            |
| <b>Immunogen</b>           | A synthetic phosphorylated peptide around S7 of human CENPA (NP_001035891.1). |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.                    |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                                   |
| <b>Note</b>                | For research use only                                                         |
| <b>Isotype</b>             | IgG                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                    |
| <b>Purity</b>              | Affinity purification                                                         |
| <b>Uniprot ID</b>          | <b>P49450</b>                                                                 |
| <b>Entrez</b>              | <b>1058</b>                                                                   |
| <b>Dilution Range</b>      | WB 1:500 - 1:2000                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                               |