

Product Datasheet

EGFR (Phospho-Y1173) antibody (orb1151137)

Description	Rabbit polyclonal antibody to EGFR (Phospho-Y1173)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y1173 of humanEGFR
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	175KDa
Uniprot ID	P00533
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.