

Product Datasheet

GLIPR1L2 antibody (orb1151143)

Description	Rabbit polyclonal antibody to GLIPR1L2
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GLIPR1L2(NP_689649.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q4G1C9
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.