

Product Datasheet

mTOR (Phospho-S2448) antibody (orb1151164)

Description	Rabbit polyclonal antibody to mTOR (Phospho-S2448)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic phosphorylated peptide around S2448 of human mTOR (NP_004949.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	289KDa
Uniprot ID	P42345
Entrez	2475
Dilution Range	WB 1:500 - 1:2000 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.