

Product Datasheet

RAP1GDS1 antibody (orb1151191)

Description	Rabbit polyclonal antibody to RAP1GDS1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RAP1GDS1 (NP_066982.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	66kDa
Uniprot ID	P52306
Entrez	5910
Dilution Range	WB 1:200 - 1:2000
Expiration Date	12 months from date of receipt.