

Product Datasheet

TBC1D16 antibody (orb1151211)

Description	Rabbit polyclonal antibody to TBC1D16
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Recombinant fusion protein of human TBC1D16(NP_061893.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TBP0
Dilution Range	IHC 1:50 - 1:100IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.