

Product Datasheet

HOMER1 antibody (orb1151217)

Description	Rabbit polyclonal antibody to HOMER1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic Peptide of human HOMER1
Preservatives	Buffer: PBS with 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q86YM7
Entrez	9456
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.