

Product Datasheet

[KO Validated] MCL1 antibody [Out of stock] (orb1151224)

Description	Rabbit polyclonal antibody to [KO Validated] MCL1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human MCL1(NP_068779.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35KDa
Uniprot ID	Q07820
Dilution Range	WB 1:500 - 1:1000 IHC 1:50 - 1:100 IF 1:50 - 1:100
Expiration Date	12 months from date of receipt.