

Product Datasheet

Tau (Phospho-T205) antibody (orb1151225)

Description	Rabbit polyclonal antibody to Tau (Phospho-T205)
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T205 of human Tau(NP_001190181.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	50-80KDa
Uniprot ID	P10636
Dilution Range	WB 1:500 - 1:2000 IHC 1:50 - 1:100
Expiration Date	12 months from date of receipt.