

Product Datasheet

[KO Validated] LC3B antibody [Out of stock] (orb1151289)

Description	Rabbit polyclonal antibody to [KO Validated] LC3B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human LC3B(NP_073729.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa/16kDa
Uniprot ID	Q9GZQ8
Dilution Range	WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.