

Product Datasheet

NeuN antibody (orb1151303)

Description	Rabbit polyclonal antibody to NeuN
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human NeuN(NP_001076045.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46-55KDa
Uniprot ID	A6NFN3
Dilution Range	WB 1:500 - 1:2000IHC 1:50 - 1:200
Expiration Date	12 months from date of receipt.