

Product Datasheet

STING (Phospho-S365) antibody (orb1151318)

Description	Rabbit polyclonal antibody to STING (Phospho-S365)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S365 of mouse STING.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	41KDa
Uniprot ID	Q3TBT3
Entrez	72512
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.