

Product Datasheet

NGF antibody (orb1151335)

Description	Rabbit polyclonal antibody to NGF
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NGF(NP_002497.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27KDa
Uniprot ID	P01138
Dilution Range	WB 1:500 - 1:1000
Expiration Date	12 months from date of receipt.