

Product Datasheet

SARS-CoV-2 Spike antibody (orb1151340)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Spike
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	DOT, ICC, IF, IP, WB
Immunogen	A synthetic peptide of coronavirus Spike(YP_009724390.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110KDa/180KDa
Uniprot ID	P0DTC2
Dilution Range	WB 1:500 - 1:2000IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.