

Product Datasheet

SARS-CoV-2 N Protein antibody (orb1151360)

Description	Rabbit polyclonal antibody to SARS-CoV-2 N Protein
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of coronavirus N Protein(NP_828858.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55KDa
Uniprot ID	P59595
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.