

Product Datasheet

Acetyl-Histone H4-K16 antibody (orb1151381)

Description	Rabbit polyclonal antibody to Acetyl-Histone H4-K16
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic acetylated peptide around K16 of human Histone H4(NP_003539.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	11kDa
Uniprot ID	P62805
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.