

Product Datasheet

Ubiquityl-Histone H2A-K119 antibody (orb1151426)

Description	Rabbit polyclonal antibody to Ubiquityl-Histone H2A-K119
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human HIST1H2AB.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	17KDa
Uniprot ID	P04908
Entrez	30128335
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.