

Product Datasheet

p120 Catenin antibody (orb1151677)

Description	Rabbit monoclonal antibody to p120 Catenin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A portion of amino acids 900-1,000 was used as the immunogen for the recombinant Delta 1 Catenin antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the recombinant Delta 1 Catenin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant Delta 1 Catenin antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	O60716
Hazard Information	This recombinant Delta 1 Catenin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 2-4ug/ml
Expiration Date	12 months from date of receipt.