

Product Datasheet

ACE antibody (orb1151801)

Description	Mouse monoclonal antibody to ACE
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full-length recombinant human ACE/CD143 protein was used as the immunogen for the Angiotensin Converting Enzyme antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the Angiotensin Converting Enzyme antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the Angiotensin Converting Enzyme antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein A/G affinity
Uniprot ID	P12821
Hazard Information	This Angiotensin Converting Enzyme antibody is available for research use only.
Dilution Range	ELISA (order BSA-free format for coating)
Expiration Date	12 months from date of receipt.