


Product Datasheet

Phospho-PLC gamma 1 (Tyr775) antibody (orb1151958)

Description	Phospho-PLC gamma 1 (Tyr775) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 1 around the phosphorylation site of Tyr775 (AL(p-Y)EG)
Target	Phospho-PLCG1 (Tyr775)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	148 kDa
Uniprot ID	P19174
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat brain
using
Phosph...