



Product Datasheet

CD21 antibody (orb1151979)

Description	CD21 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD21
Target	CD21
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	111 kDa
Uniprot ID	P20023
Dilution Range	WB=1:500-2000, IHC-P=1:1000, ICC=1:100
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human
spleen using
Pri...



Western
blot:Primary
antibody at
1:2000 ...