

Product Datasheet

LYPD3 antibody (orb1151980)

Description	LYPD3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human LYPD3
Target	LYPD3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	33 kDa
Uniprot ID	O95274
Dilution Range	WB=1:500-1000, IHC-P=1:1000, ICC=1:50, IF=1:50-100
Expiration Date	12 months from date of receipt.