



Product Datasheet

WNT10B antibody (orb1151991)

Description	WNT10B antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human WNT10B
Target	WNT10B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1/Kappa
Clonality	Recombinant
MW	40 kDa
Uniprot ID	O00744
Dilution Range	WB=1:500-2000, IHC-P=1:50-1:100, ICC=1:50
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human liver
cancer usi...



IHC-P
analysis of
rat brain
using
WNT10B...