



Product Datasheet

alpha Actinin 4 antibody (orb1152016)

Description	alpha Actinin 4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human alpha Actinin 4
Target	alpha Actinin 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	105 kDa
Uniprot ID	Q43707
Dilution Range	WB=1:200-1000, IHC-P=1:20-1:100, ICC=1:20-1:100
Expiration Date	12 months from date of receipt.



Western blot: Primary antibody at 1:5000 ...



Western blot: Lane 1: Mouse Skin; Lane 2: Mo...