

Product Datasheet

Human Tau-4 antibody (orb1152038)

Description	Human Tau-4 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Tau-4 protein (1-441/441 aa)
Target	Human Tau-4
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46/79 kDa
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.