

## Product Datasheet

### Human Tau-4 antibody (orb1152039)

|                            |                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Tau-4 antibody                                                                |
| <b>Species/Host</b>        | Mouse                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                        |
| <b>Tested Applications</b> | ELISA, WB                                                                           |
| <b>Immunogen</b>           | Recombinant Human Tau-4 protein (1-441/441 aa)                                      |
| <b>Target</b>              | Human Tau-4                                                                         |
| <b>Preservatives</b>       | 0.01M PBS(pH7.4)                                                                    |
| <b>Form/Appearance</b>     | Liquid                                                                              |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year.<br>Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                               |
| <b>Isotype</b>             | IgG1                                                                                |
| <b>Clonality</b>           | Monoclonal                                                                          |
| <b>MW</b>                  | 46/79 kDa                                                                           |
| <b>Uniprot ID</b>          | <b>P10636</b>                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                     |