



Product Datasheet

Prostate Specific Antigen antibody (orb1152046)

Description	Prostate Specific Antigen antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Prostate Specific Antigen (PSA) Antigen, Native
Target	Prostate Specific Antigen
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	23 kDa
Uniprot ID	P07288
Dilution Range	WB=1:500-2000, IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
prostate
using P...



Western
blot: Lane
1: Prostate
Specific
An...