

Product Datasheet

phospho-Histone H3 (Ser10) antibody (orb1152070)

Description

phospho-Histone H3 (Ser10) antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human Histone H3 around the phosphorylation site of Ser10 (RK(p-S)TG)

Target

phospho-Histone H3 (Ser10)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

15 kDa

CAS Number

Histone H3.1

Uniprot ID
P68431
Dilution Range

WB=1:500-2000

Expiration Date

12 months from date of receipt.