


Product Datasheet

CH25H antibody (orb1152120)

Description	CH25H antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human CH25H protein
Target	CH25H
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	O95992
Dilution Range	IHC-P=1:400-800, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat brain
using CH25H
...