

Product Datasheet

RPE65 antibody (orb1152131)

Description	RPE65 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse RPE65
Target	RPE65
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	Q91ZQ5
Dilution Range	IHC-P=1:1000-4000, WB=1:500-2000
Expiration Date	12 months from date of receipt.