


Product Datasheet

TurboGFP antibody (orb1152139)

Description	TurboGFP antibody
Species/Host	Mouse
Reactivity	Tag
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant TurboGFP protein expressed in HEK293T cell.
Target	TurboGFP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	85 kDa
Dilution Range	WB=1:50000-200000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1:2: Recombinant Tur...