

Product Datasheet

phospho-Tau (Thr217) antibody (orb1152152)

Description	phospho-Tau (Thr217) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr217 (LP(p-T)PP)
Target	phospho-Tau (Thr217)
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	79 kDa
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.