

Product Datasheet

**SARS-CoV-2 (2019-nCoV) Nucleocapsid protein,
Biotin conjugated antibody (orb1152154)**

Description

SARS-CoV-2 (2019-nCoV)

	conjugated antibody
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protei
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	≥1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	P0DTC2
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558