

## Product Datasheet

### RING2/RING1B/RNF2 Antibody (orb1152320)

## Description

RING2/RING1B/RNF2 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, IHC, WB

### Applications

### Immunogen

E.coli-derived human RING2/RING1B/RNF2 recombinant protein (Position: M1-K336).

### Form/Appearance

Lyophilized

### Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

### Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

### Note

For research use only

### Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

### Isotype

Rabbit IgG

### Clonality

Polyclonal

### MW

38 kDa

### Uniprot ID


[Q99496](#)

### Dilution Range


Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Direct ELISA, 0.1-0.5 µg/ml, Human

### Expiration Date


12 months from date of receipt.



IHC analysis of RING2/RING1B/RNF2 using ...



IHC analysis of RING2/RING1B/RNF2 using ...



IHC analysis of RING2/RING1B/RNF2 using ...