


Product Datasheet

Integrin beta 3/ITGB3 Antibody (orb1152336)

Description

Integrin beta 3/ITGB3 Antibody



Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IHC, WB

Applications

Immunogen

E.coli-derived human Integrin beta 3/ITGB3 recombinant protein (Position: N29-Q523).

IHC analysis
of Integrin
beta 3/ITGB3
us...

Form/Appearance


Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.



IHC analysis
of Integrin
beta 3/ITGB3
us...

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


100 kDa

Uniprot ID

[P05106](#)

Dilution Range

Western blot, 0.25-0.5 µg/ml, Human, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Rat Flow Cytometry, 1-3 µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml, Human



IHC analysis
of Integrin
beta 3/ITGB3
us...

Expiration Date

12 months from date of receipt.