

## Product Datasheet

### HLA-DR/HLA-DRA Antibody (monoclonal, 5B13F7) (orb1152388)

## Description

HLA-DR/HLA-DRA Antibody (monoclonal, 5B13F7)

## Species/Host

Mouse

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

FC, IHC, WB

## Applications

## Immunogen

E.coli-derived human HLA-DR/HLA-DRA recombinant protein (Position: I26-L254).

## Form/Appearance

Lyophilized

## Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

## Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

## Note

For research use only

## Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

## Isotype

Mouse IgG2b

## Clonality

Monoclonal

## MW

35-37 kDa

## Uniprot ID



**P01903**

## Dilution Range


Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10<sup>6</sup> cells, Human

## Expiration Date


12 months from date of receipt.



IHC analysis of HLA-DRA using anti-HLA-D...



IHC analysis of HLA-DRA using anti-HLA-D...



IHC analysis of HLA-DRA using anti-HLA-D...