

Product Datasheet

Drebrin/DBN1 Antibody (monoclonal, 4F6E7) (orb1152404)

Description

Drebrin/DBN1 Antibody (monoclonal, 4F6E7)

Species/Host

Mouse

Reactivity

Human, Monkey, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IHC, WB

Applications

Immunogen

E.coli-derived human Drebrin/DBN1 recombinant protein (Position: H9-D649).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2a

Clonality

Monoclonal

MW

120 kDa

Uniprot ID


[Q16643](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat, Monkey Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Rat Flow Cytometry, 1-3 µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



IHC analysis of DBN1 using anti-DBN1 ant...



IHC analysis of DBN1 using anti-DBN1 ant...



IHC analysis of DBN1 using anti-DBN1 ant...