

Product Datasheet

NDFIP2 Antibody (monoclonal, 10D6D7) (orb1152407)

Description

NDFIP2 Antibody (monoclonal, 10D6D7)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications

Immunogen

E.coli-derived human NDFIP2 recombinant protein
(Position: M16-L336).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration
of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After
reconstitution, at 4°C for one month. It can also be
aliquotted and stored frozen at -20°C for six months.
Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive
results. Other applications have not been tested. Optimal
dilutions should be determined by end users. Adding 0.2
ml of distilled water will yield a concentration of 500
µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

39 kDa

Uniprot ID


[Q9NV92](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human
Immunohistochemistry(Paraffin-embedded Section), 2-5
µg/ml, Human
Immunocytochemistry/Immunofluorescence, 5 µg/ml,
Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



IHC analysis
of NDFIP2
using anti-
NDFIP2...



IHC analysis
of NDFIP2
using anti-
NDFIP2...



IHC analysis
of NDFIP2
using anti-
NDFIP2...