


Product Datasheet

Mouse IL-27 p28 ELISA Kit (DIY Antibody Pairs) (orb1152409)

Description	Mouse IL-27 p28 ELISA Kit (DIY Antibody Pairs)
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F29-S234
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q8K316
Expiration Date	Please enquire.



Mouse IL-27 p28 EZ Set ELISA Kit standar...