


Product Datasheet

**Biotinylated Human IL-18 R1 / CD218a Protein (SPR
& MALS verified) (orb1152438)**

Description

Biotinylated Human IL-18 R1, Fc, is


Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

IL-18 R1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc & C-Avi

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 63.9 kDa. The protein migrates as 70-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Biotinylated Human IL-18 R1, Fc, (orb1152438) is expressed from human 293 cells (HEK293). It contains AA Ala 19 - Arg 329 (Accession # AAH69575.1).

MW

63.9 kDa

NCBI

[AAH69575.1](#)

Expiration Date

6 months from date of receipt.