



Product Datasheet

TSARG4 antibody (orb11525)


Description	TSARG4 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TSARG4 (161-275/379 aa)
Target	SPAG4L/SRG4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	Q8TC36
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



IHC-P of rat pancreas tissue (dilution a...



Western blot analysis of rat brain lysat...