

Product Datasheet

CD41 antibody (orb1152520)

Description	Mouse monoclonal antibody to CD41
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	CD41 is the integrin alpha chain 2b forming a fibrinogen receptor in heterodimeric complex with beta 3.
Target	CD41
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P08514
Expiration Date	12 months from date of receipt.