

Product Datasheet

Rat monoclonal anti-mouse Interleukin-6 (IL-6) (orb1152561)

Description	Rat monoclonal anti-mouse Interleukin-6 (IL-6)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP
Immunogen	Mouse recombinant IL-6
Storage	Reconstitute the antibody with sterile PBS and the reconstituted antibody can be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles. Lyophilized samples are stable for 2 years from date of receipt when stored at -70°C.
Note	For research use only
Isotype	Rat IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.