

Product Datasheet

Rat Monoclonal anti-mouse Progranulin (orb1152644)

Description	Rat Monoclonal anti-mouse Progranulin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant Mouse Progranulin
Storage	Reconstitute the antibody with 500ul sterile PBS and the final concentration is 200ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2
Clonality	Monoclonal
Dilution Range	WB (1 to 1000-2000), ICC (1 to 100- 200)
Expiration Date	12 months from date of receipt.