

Product Datasheet

Rat Monoclonal anti-mouse IL-6 Receptor (orb1152921)

Description	Rat Monoclonal anti-mouse IL-6 Receptor
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Mouse recombinant protein of the extracellular domain of IL-6 R
Storage	Reconstitute the antibody with 200 ul sterile PBS (final concentration 500ug/ml). Reconstituted antibody can be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles
Note	For research use only
Isotype	Rat IgG2
Clonality	Monoclonal
Dilution Range	WB (1 to 500-1000), IHC-Fr (1 to 50-100)
Expiration Date	12 months from date of receipt.