

Product Datasheet

Rat Monoclonal anti-mouse CD45 (orb1152975)

Description	Rat Monoclonal anti-mouse CD45
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	Cell lysates of mouse leukocytes
Storage	Reconstitute the antibody with 500 ul sterile PBS (final concentration 400ug/ml). Reconstituted antibody can be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2
Clonality	Monoclonal
Dilution Range	IHC-P (1 to 100-400), IP (1 to 100)
Expiration Date	12 months from date of receipt.