

Product Datasheet

Rat Monoclonal anti-mouse CD86 (orb1152989)

Description	Rat Monoclonal anti-mouse CD86
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	LPS stimulated Mouse B-Cells
Storage	Reconstitute the antibody with 400 ul sterile PBS (final concentration 500ug/ml). Reconstituted antibody can be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2
Clonality	Monoclonal
Dilution Range	IHC-P (1 to 300-2000), FC (1 to 100-1000)
Expiration Date	12 months from date of receipt.